

Options Tricycles



Velcro footplates
(small model)



crutch holder



T-shaped saddle pin



Velcro pedals



rectangular
steering wheel
(large model)



kilometer counter



drum brake, built into
the front wheel axle



front light



accompanying bar



footclasp



accompanying bar with
hip support and belt



rear light

Technical information Tricycles

	length	width	height saddle		height steering wheel		diameter
			Min.	Max.	Min.	Max.	
2102	110	67	48	68	72	84	30
2103	120	67	54	74	81	96	40
2115	135	66	56	70	75	85	40
2117	162	67	75	90	96	100	50
2119	180	80	85	105	95	110	60

Sizes in cm

Deviations in sizes (1,5 cm) are possible

Your guarantee of quality

Regarding quality, the Express 2000 complies with the **highest norms laid down by certified test institutes** in most European countries.

Adaptable to your needs

Vermeiren's range of quality tricycles is the result of continuous research and more than 45 years of experience. Each model is thoroughly tested and evaluated. One can choose from a wide range of options. Your tricycle must be adapted to you, not the other way around.

N.V. Vermeiren N.V.


Vermeirenplein 1-15
B-2920 Kalmthout, Belgium

Tel.: +32 (0)3 620 20 20
Fax: +32 (0)3 666 48 94
Visit our website: www.vermeiren.com



VERMEIREN
Opening up new horizons in mobility

For more information,
contact your local dealer.



Subject to technical alterations

Tricycles

Mobility with total freedom



2119



Tricycles

Mobility is essential

Freedom of movement is important for everyone and almost equals independence. Furthermore, mobility and even more so cycling are important parts of rehabilitation.

That's why Vermeiren offers a wide range of tricycles. Five different models, each with its own characteristics and features, the result of more than 45 years of research and experience.



2102



2117

Tricycles, their characteristics and possible options

The table below gives you a full overview of the various models, including detailed specifications and information about the available colours and options.



footplate with legclamps



adjustable pedal bars



broad hip support with belt

	01	02	04	05	06	08	09	10	11	12	14	15	16	17	18	19	20	21	22	23	24	25	colours	
2102	<input type="checkbox"/>	<input type="checkbox"/>			s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>				<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>	red/black	
2103	<input type="checkbox"/>	<input type="checkbox"/>			s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>				<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>	yellow/black	
2115	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	s	s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>	yellow/black
2117	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	s	s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					blue
2119	<input type="checkbox"/>	<input type="checkbox"/>	s	s	s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	s		grey

s = standard = option

- 01 neutral
- 02 back pedal brake
- 04 front and rear light
- 05 drum brake, built into the front wheel axle
- 06 fixed sprocket
- 08 T-shaped saddle pin allows the saddle to be adjusted back and forward
- 09 Velcro pedals
- 10 footclasp
- 11 Velcro footplates (small model)
- 12 Velcro footplates (large model)
- 14 broad hip support with belt

- 15 crutch holder
- 16 rectangular steering wheel (small model)
- 17 rectangular steering wheel (large model)
- 18 kilometer counter
- 19 adjustable pedal bars, in order to position the pedals higher
- 20 footplate with legclamps
- 21 parking brake
- 22 back pedal brake, 3-speed
- 23 differential
- 24 accompanying bar
- 25 accompanying bar and hip support



2115



2103